


Enclosed are your child's newly created earmold(s). The length of the tubing can be easily trimmed to a perfect fit with these Instructions....



1. First, place the earmold in your child's ear and ensure a proper fit. Then, hang the hearing aid over the back of the ear.
2. Next, measure with a pen or marker and mark SLIGHTLY ABOVE the point at which the earmold tubing and the earhook meet.
3. Finally, remove the earmold and aid from the child's ear, and make a blunt, straight cut ABOVE the marked line- Always cut above the line, so that cutting the tubing too short won't occur.



TIPS:

*If the aid is sitting too far above the ear when attached to the earmold, trim the top of the tubing, just a little at a time, until you reach the optimal fit.

*The tubing may be difficult to apply to the earhook of the aid at first; gently wiggling the tubing onto the earhook will typically work .

*Don't worry - if too much tubing is removed when trimming, a new earmold is not necessary; we can easily re-tube the earmold in most cases.

Please call or e-mail if you have any questions or concerns.

Thank You!